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fNIRS / FUNCTIONAL NEAR-INFRARED SPECTROSCOPE

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This equipment (Fig.1) allows recording the functional activity of the brain by optical imaging technique. It enables to study the spatial and temporal brain activations for basic research in the field of vision. This technique is a type of functional neuroimaging technology that offers a non-invasive method of monitoring of brain activity (Fig.3).

Figure 1
Experience lors de la session de formation le 23 avril 2013 à L'EquipEx IrDIVE
Experiment during the training the 23th of May 2013 at the EquipEx IrDIVE

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PRINCIPE

The principle of fNIRS is to measure changes of absorption in near-infrared light between oxy-hemoglobin (Oxy-Hb) and deoxy-hemoglobin (Deoxy-Hb) in blood flow on the surface of brain (Fig.3). An activation of brain increases blood flow and vascular expansion, as a result Oxy-Hb increases when the desoxy-Hb decreases. 3 laser sources are used in the fNIRS: 780nm, 805 and 830nm, they are transmitted by optical fibres (*transmitter fibres*) set on the head, then absorbed in part by the tissues. The transmitted light signal is collected by other optical fibres (*receiver fibres*), and then the signal is analysed and processed (Fig.4) to represent the absorption of *oxy-Hb* and *desoxy-Hb*. Different parts of brain are activated; it is then possible to make a map of brain.


Figure 2
Emission et réception de la lumière.
Emission and reception of light (from Shimadzu training).

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SPECIFICATIONS

- 48 mesures points
- Wavelengths of laser sources: 780nm, 805 and 830nm
- Synchronization between Command computer and recording computer
 - Real time measures
 - Able to work with 64 channels EEG
 - Real time Display of active areas
 - Limited Exterior electromagnetic perturbations.

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APPLICATIONS

- Studies of neural basis of perception and cognition.
- Studies of brain areas involved in decision and organization of voluntary acts.
- Studies of the correlation between the hemodynamic activities (fNIRS recorded) and electroencephalographic activities (EEG recorded) in a motor planning task.

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